

J.K.

6350L Fluorescent Immunohistochemistry Protocol 9/17/22

20% PFA (50mL)

$$20\% = \frac{20g}{100mL} = \frac{xg}{50mL} \Rightarrow \boxed{x = 10g}$$

4% PFA (50mL)

$$20\% \times mL = 4\% 50mL \Rightarrow \boxed{x = 10mL}$$

20% Tween-20 (50mL)

- Measure <50mL of 1x PBS into 50mL conical tube
- Add 10g (if solid) or (10mL) (if liquid) of Tween-20.

Note! - From 20% calculation above

PBST (1L) $\checkmark 100$

- Measure ~~500~~ 500mL of 10x PBS into 500mL G.C.
- Add 1mL of 20% tween-20(aq.)

↳ 0.1% of PBST is tween-20.

$$\hookrightarrow 0.001 \times 1000mL = 1mL$$

2.5 M Glycine stock (50 mL)

- Can measure $> 50 \text{ mL}$ H_2O to dissolve glycine in.
- $2.5 \text{ M} = \frac{2.5 \text{ mol}}{1 \text{ L}} = \frac{x \text{ mol}}{.050 \text{ L}} \Rightarrow x \text{ mol} = 0.125 \text{ mol}$
- $0.125 \text{ mol} \times 75.07 \text{ g/mol} = \boxed{9.38 \text{ g}}$ \rightarrow mw of glycine
- Bring to Volume: 50 mL .

100 mM Sodium Citrate Stock 500 mL

- $100 \text{ mM} = \frac{0.1 \text{ mol}}{1 \text{ L}} = \frac{x \text{ mol}}{0.500 \text{ L}} \Rightarrow x \text{ mol} = 0.05 \text{ mol}$
- $0.05 \text{ mol} \times 258.06 \text{ g/mol} = \boxed{12.9 \text{ g}}$